

**BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAMME
UNDER
NATIONAL HEALTH MISSION**



CYRIX
HEALTHCARE PVT LTD

No. : 150954

SERVICE PROVIDER
Tender No. WO-37/2021-2022/698

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

Service Report

Health Facility ... G.H. Kasaragod ...

Address ... Kasaragod ...

Ph : 98466 22912

Call Registration Date : 27-10-23
 Caller ID : 77032
 Date of Visit : 27-10-23
 Asset No. : 1410927
 EQPT Name : debi baillator
 Manufacture Philips Model : Heart start xl
 S. No. U.S.O.0578.196 Dept. ICU

Service Classification : Breakdown Call PMS Calibration Cust. Training

Problem Identified : Not switching on

Action Taken : checked the debi baillator and found problem with the snps pcb and main pcb, machine installed on 9-1-10 and aged up to 13+ years. OEM declared that the Philips debi baillator model Heart start xl is an end of support device hence spare and service are not available. Some recommendation given for condensation.

Completed Date : 27-10-23 Time : 4:00 Spare Required

Description	Qty.	Part Number	PR Number		
			Date	Start Time	End Time
1.					
2.					
3.					
Cyrix Engineer			27-10-23	3:00	4:00
Ajith Babu K					

Customer Remark Completed Pending

Service Engineer Name : Ajith Babu K
Signature : [Signature]
Date : 27-10-23
Contact Number : 8848079288


Customer Name : JEEPH - C
Signature : [Signature]
Date : 23.10.23
Contact Number : 9469290216
Designation : [Signature]
Hospital Seal : [Seal]
 - 671 1

**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAM (BEMP)**

Recommendations for Beyond Economic Repair (BER)

PROFORMA

Sl. No	Particulars	Details
1	Name of District	Kasaragod
2	Name of Hospital	SH Kasaragod
3	Name of Equipment with Make, Model and Serial Number	Defibrillator Philips, Heartstart XL U500578196
4	Equipment ID & Barcode	1410927
5	Date of purchase / Year of manufacture / Installation Date	9-1-2010
6	Warranty details (Yes/No)	NO
7	*AMC/ CAMC Period agreed at the time of purchase	NO
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	27-10-23
9	Action taken	Attached cyrize service report
10	Present status of the equipment (Fully damaged / partially damaged)	fully damaged
11	Recommendations for repair (required service details)	not recommending job repair
12	Cost of spares (specify parts and cost)	NA

13	Asset Value	177,000/-
14	* Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	NA
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	Attached cyrix service report
16	Reasons for recommending the equipment as BER	checked the machine and found problem with the main PCB and SMPS PCB
17	Name & Signature of CYRIX Authority	Ajith Babu. K 

*Not mandatory #Based on the period of life and value as per the BER guidelines

* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:

Checked the machine and found problem with SMPS and main PCB. Defibrillator installed on q1110. OEM declared that, defibrillator model Heartstart XL as an end of support device. OEM report attached.

Rev. 31/10/23
 Signature of JC BM (NHM) **Jr. CONSULTANT**
 Biomedical Engineering
 Areyvakeralam (NHM)
 Kozhikode, Kozhikode

Date



Signature of
 Superintendent / Medical Officer (i/c)


 സൂപ്രണ്ട്

ഗവ. ജനറൽ ആശുപത്രി
 കാസർഗോഡ്
 കാസർഗോഡ് പി.ഒ. - 671 121

1st Floor, "Vrindavan", Door No. 43/1606 B,
North Railway Station Road, Ernakulam, Cochin - 682 018, India


Dare: 6th Apr 2022

To Whomsoever it may concern

This is to inform that Philips India Ltd have declared defibrillator model Heartstart XL as an end of support device. Hence spares and service are not available for this model.

Thanking you,

Yours Faithfully


Davis George
Manager
Anamdev Engineers



Costs cost ✓

Cardiac Defibrillator & Monitor DF 2509, IR

B/f from old book p. No - 31

27/4/09 condemnation Bd from store 1 - 1 ✓
Rs 11000/- (BPL)

2/4/07 Londonmachin Bd from store 1 - 2 ✓
Model - BPL

DF - 2509 R Rs 1,80,000/-
D.No - 205 (IRN)

10/7/08 Bd - from store 1 - 3 ✓
Silicon lab pvt ltd - G. Care
Chennai cost 1.


Ph: - 91-44-22394566
91-44-22394652
91-44-22396762

(Warranty card not Bd)

9/1/10 Bd from store 1 - 4 ✓
S.No - 00578196
Rs - 1,84,151/-

Accessories in store - Bd
Philips - Hertostat
cost not mentioned.

6/12/12 Transferred to unseizable register 2 - 1 ✓
page no. 32.
condemnation done on 4/3/21





REF 861266
 SN US00578196
 OPT RLJ13

LCD DISPLAY
 PMS M735R


Service # 861266
 SN US00578196



100 - 240 V
 50 - 60 Hz
 1.5 A



Production Order
 2009-12



注册证号: SFDA(1)20083212524
 生产企业: Manufactured by

产品标准编号: YZB/USA4546-2006

注册地址 / 生产地址: Philips Medical Systems
 3000 Minuteman Road
 Andover, MA 01810 USA 美国
 联系方式: Tel 8008100038 (中国)



Small text at the bottom of the label, including manufacturing details and company information.

ECG In

ECG Out

DANGER: Possible explosion hazard if used in the presence of flammable anesthetics or concentrated oxygen.

WARNING: Hazardous electrical output. This equipment is for use only by qualified medical personnel.

NOTICE: Plug into rms AC line during standby periods.

MANUAL OPERATION

DEFIBRILLATION:

- 1 Select Energy.
- 2 Press CHARGE.
- 3 Press SHOCK.



ECG MONITORING:

- Select Manual On.
- Apply Electrodes.
- Press Lead Select.

AED OPERATION

DEFIBRILLATION:

- 1 Apply Pads. Select AED On.
- 2 Follow prompts. If instructed, Press ANALYZE.
- 3 If shock advised, press SHOCK.

CYRIX BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME (BEMP)

TOLL FREE NO: 1800 - 425-7669

BAR CODE - (8004-890615225)
1410927



RBER

CYRIX

Ticket No: 33022
 RBER approved on: 27/08/21
 Condemn on:
 Remark: not working properly

Mark

