

*Handwritten:* 8-18  
*Handwritten:* KMSCL

**CYRIX** BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME (BEMP)  
TOLL FREE NO: 1800 - 425-7669  
BAR CODE - (8004-890615225) 0540460

KMSCL KTPL  
TOLL FREE NO:- 1800-270-0690  
(8004) 8907579 - 70051036



Imported & Marketed By:-  
Doctor Bellam & Sons Pvt. Ltd.  
Head Office (Customer Care): 2/17 Ann  
Customer Care Ph: 011-23270270, Custom  
Model No: Progress  
Description: Bathroom weighing scale  
Quantity: 1 Number  
Month & Year of Import: July ' 2018  
Mfg: Dec 2009. Inclusive Of All Taxes

SPECIFICAT  
Size  
Capacity  
Weight

LED Display 1.3 inches  
Automatic Zero Reset  
Low Battery Indicator



# BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME UNDER NATIONAL HEALTH MISSION



**SERVICE PROVIDER**  
Tender No. WO-37/2021-2022/698

## CYRIX

HEALTHCARE PVT LTD

**No. :** 172414

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

**Service Report**

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala  
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

<p>Health Facility : <u>PHC</u></p> <p>Address : <u>PHC Thiruvapur</u> <u>Kottayam</u></p> <p>Ph : <u>9495244316</u></p>	<p>Call Registration Date : <u>13/12/2023</u></p> <p>Caller ID : <u>85671</u></p> <p>Date of Visit : <u>16/12/2023</u></p> <p>Asset No. : <u>0540460</u></p> <p>EQPT Name : <u>Weighing Machine</u></p> <p>Manufacture <u>Sanso</u> Model : <u>ni</u></p> <p>S. No. <u>ni</u> Dept. <u>opd</u></p>
------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Service Classification :** Breakdown Call  PMS  Calibration  Cust.Training

**Problem Identified :** Weighing machine not working

**Action Taken :** Checked and found battery broke and leaked. The battery connectors, internal body parts and connectors are tested. Recommended for replacement.

Completed  Date : ..... Time : ..... Spare Required

Spare Replaced  Requested

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrux Engineer	Date	Start Time	End Time
<u>Sanish Sivan</u>	<u>16/12/23</u>	<u>10:00am</u>	<u>10:30am</u>

**Customer Remark** Completed  2 Pending

*(Handwritten signatures and notes)*


<p><b>Service Engineer Name :</b> <u>Sanish Sivan</u></p> <p><b>Signature :</b> <i>(Signature)</i></p> <p><b>Date :</b> <u>16/12/2023</u></p> <p><b>Contact Number :</b> <u>9946646571</u></p>	<p><b>Customer Name :</b> <u>Dr. RAJESH</u></p> <p><b>Signature :</b> <i>(Signature)</i></p> <p><b>Date :</b> <u>16/12/23</u></p> <p><b>Contact Number :</b> <u>9495058737</u></p> <p><b>Designation :</b> <u>Dr. Am Sivan</u></p> <p><b>Hospital Seal :</b> <i>(Seal)</i></p>
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT  
MAINTENANCE PROGRAM (BEMP)**

**Recommendations for Beyond Economic Repair (BER)**

**PROFORMA**

Sl. No	Particulars	Details
1	Name of District	Kottayam
2	Name of Hospital	PHC Thiruvapur
3	Name of Equipment with Make, Model and Serial Number	Weighing Machine Make: Samso Model: nil S/n: nil
4	Equipment ID & Barcode	0540460
5	Date of purchase / Year of manufacture / Installation Date	01-10-2018
6	Warranty details (Yes/No)	Nil
7	*AMC/ CAMC Period agreed at the time of purchase	Nil
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	13-12-2023
9	Action taken	checked the equipment and found not working
10	Present status of the equipment (Fully damaged / partially damaged)	Fully damaged.
11	Recommendations for repair (required service details)	checked and found battery broke and leaked. Connectors are rusted and broken. Internal connectors and body parts are wear and tear. Repairing / replacement of battery connectors and internal parts.
12	Cost of spares (specify parts and cost)	Spares not available


13	Asset Value	2590/-
14	*Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	Spare not available
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	cyrix report recommends the equipment to be replace.
16	Reasons for recommending the equipment as BER	The repairing of the equipment is economically not feasible. Hence recommended the equipment as BER
17	Name & Signature of CYRIX Authority	 Sanis Swas

\*Not mandatory  #Based on the period of life and value as per the BER guidelines


\* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:

Asset Value - 2590  
 Cost of spare = Spare not available  
 Installation Date = 11/10/2018  
 Life period = 5 year 2 months.  
 As per the BER guidelines recommended as BER

  
 CONSULTANT BIOMEDICAL ENGINEER  
 NATIONAL HEALTH MISSION  
 KOTLAYAM  
 Signature of JC BM (NHM)

Date

  
 Signature of Superintendent / Medical Officer (i/c)



11.15.20  
 Dr. N. N. N.  
 MEDICAL OFFICER  
 PHC Thiruvappu  
 Kotlayam



SAMSO

SPECIFICATIONS:  
Item Code:  
Weight:  
Metric Number:  
© S. M. 2001



WRIGHTWORLD MACHINE  
(DIGITAL)

<p>1/10/18 (Ss Andacio Page:18)</p>	<p>Received from karunya medicine depot as per 1000 11457 dt 1/10/18.</p>	<p>2</p>	<p>2</p>	<p>2</p>	<p>2590/- 1</p>
<p>1/10/18</p>	<p>issued to s/n</p>	<p></p>	<p>2</p>	<p>NIL Rd</p>	<p>9.00</p>