


**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAM (BEMP)**

Recommendations for Beyond Economic Repair (BER)

PROFORMA

Sl. No	Particulars	Details
1	Name of District	Kottayam
2	Name of Hospital	PHC Poonjar
3	Name of Equipment with Make, Model and Serial Number	BP Apparatus Make: Anita Industries Model: Elkonster Sn: 449392
4	Equipment ID & Barcode	0541913
5	Date of purchase / Year of manufacture / Installation Date	27-11-2018
6	Warranty details (Yes/No)	No
7	*AMC/ CAMC Period agreed at the time of purchase	No
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	25-01-2004, 92766
9	Action taken	Serviced.
10	Present status of the equipment (Fully damaged / partially damaged)	Fully damaged.
11	Recommendations for repair (required service details)	Required spare parts of filter, mercury, BP cuff, bladder and bulb with valve. Body is damaged conditions.
12	Cost of spares (specify parts and cost)	Not available


13	Asset Value	1303
14	# Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	Not available.
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	Cyrix.
16	Reasons for recommending the equipment as BER	Spare cost is almost equal to asset value. Machine is aged upto 5 years and damaged conditions. Economically not repairable.
17	Name & Signature of CYRIX Authority	Nidhish Roy 

*Not mandatory #Based on the period of life and value as per the BER guidelines

* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:


Asset value - 1323
 Cost of spare - Not available
 Life period - 5 year 2 months
 Cyrix reported that the spare not available

CONSULTANT BIOMEDICAL ENGINEER
 NATIONAL HEALTH MISSION

 Signature of JC BM (NHM)

Meningeal BP apparatus is not working

ASITHA.V.S
 SENIOR NURSING OFFICER
 CWC POONJAR

Date 30/1/2024


 Signature of Superintendent / Medical Officer (w/c)



എമ്മിഎൽ ആഫീസർ ഇൻ-ചാർജ്
 പ്രാഥമിക ആരോഗ്യകേന്ദ്രം
 പുത്തൂർ തെക്കേക്കര
 കോട്ടയം, പിൻ: 686 582



**BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAMME
UNDER
NATIONAL HEALTH MISSION**



SERVICE PROVIDER
Tender No. WO-37/2021-2022/698

CYRIX
HEALTHCARE PVT LTD

No. : 174237

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

Service Report

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

Health Facility <u>PHC Poonjar</u>	Call Registration Date : <u>25-01-2024</u>
Address <u>Poonjar</u>	Caller ID : <u>98746</u>
<u>Kottayam</u>	Date of Visit : <u>27-01-2024</u>
Ph : <u>9495188532</u>	Asset No. : <u>0541913</u>
	EQPT Name : <u>BP Apparatus</u>
	Manufacture <u>Anita</u> Model : <u>Ellcometer</u>
	S. No. <u>449398</u> Dept. <u>OP</u>

Service Classification : Breakdown Call PMS Calibration Cust.Training

Problem Identified : Low mercury, BP bladder bladder leak, cuff damaged, Required BP bulb with valve.

Action Taken : Checked. Required Filter replacement. Required mercury, BP bladder, cuff, bulb with valve. Body totally damaged condition. Machine is very old and economically not repairable.

Completed Date : 30.01.2024 Time : 10:35am Spare Required

Spare Replaced Requested

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrux Engineer	Date	Start Time	End Time
<u>Nidhish Roy</u>	<u>30-01-2024</u>	<u>10.00am</u>	<u>10.35am</u>



Customer Remark BP apparatus (mercury) not working. Completed Pending
ASITHA-V-S Senior Nursing officer.

Service Engineer Name : Nidhish Roy
Signature : [Signature]
Date : 30-01-2024
Contact Number : 9872575296

Customer Name : [Signature]
Signature : [Signature]
Date : [Signature]
Contact Number : [Signature]
Designation : [Signature]
Hospital Seal : [Signature]
കോട്ടയം, പിൻ: 686 582

BP APPARATUS

67

DATE		Rd	Id	Balance	By
	BIF from instrument - Furniture Book P. 176				
15/8/19	Installation date	Rd	2	2	
15/8/19	medicine Atomach	Put on use (ded)	2		
3/9/19	Installation date	Rd	1	3	
3/9/19	Put on use (ded)		1		
31/3/22	Installation date	Rd	1	4	
	BIF LP P. 310				
31/3/22	Put on use		1		
27/11/18	Installation date	BIF from instrument - Furniture Book			
27/11/18	P. 166	Rd	3	7	
27/11/18	Put on use		3		
17/3/10	Installation date	BIF from instrument - Furniture Book			
17/3/10	P. 146	Rd	①	8	
17/3/10	Automatic BP monitor	Digital			
	Put on use		1		
15/11/23	Installation date	Rd from Store			
	electronic LP	(can use)	1	9	
18/1/23	Put on use		1		
25/10/22	Installation date	BIF from instrument - Furniture Book P. 49			
	3 (15-16107)		⑦	16	
	4 (2914114)				
25/10/22	Transferred to unserviceable register P. 62				
1/3/23	Installation date	condemned 7 Nos as per order No.			
	390/22 of Mo FHC P. 62 dt 1/3/23				
	16 - 7 = 9	Balance	9		
25/12/23	Installation date	Transfer to Unserviceable Register			
	P. No: 62. Method, [9-1] = 8				
	(27/11/18 (many BP apparatus))				

