



# BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME UNDER NATIONAL HEALTH MISSION



**SERVICE PROVIDER**  
Tender No. WO-37/2021-2022/698

# CYRIX

HEALTHCARE PVT LTD

**No. :** 178906

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

**Service Report**

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala  
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

<p>Health Facility <u>DH Wadakkanchery</u></p> <p>Address <u>Thrissur</u></p> <p>Ph : <u>7560935204</u></p>	<p>Call Registration Date : <u>17/01/2024</u></p> <p>Caller ID : <u>91214</u></p> <p>Date of Visit : <u>19/01/2024</u></p> <p>Asset No. : <u>0812018</u></p> <p>EQPT Name : <u>DIALYSER REPROCESSOR UNIT</u></p> <p>Manufacture <u>CURA</u> Model : <u>DCS PLUS</u></p> <p>S. No. <u>06510218</u> Dept. <u>DIALYSIS</u></p>
-------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Service Classification : Breakdown Call**  **PMS**  **Calibration**  **Cust.Training**

**Problem Identified :** Not working

**Action Taken :** checked the Equipment found that Main board and Solenoid Valve let defective. Enquired for spare.

**Completed**  **Date :** ..... **Time :** ..... **Spare Required**

**Spare Replaced**  **Requested**

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrix Engineer	Date	Start Time	End Time

**Customer Remark**  **Completed**  **Pending**



**Service Engineer Name :** Tejas  
**Signature :** [Signature]  
**Date :** 19/01/2024  
**Contact Number :** 8848098422


**Customer Name :** [Signature]  
**Signature :** [Signature]  
**Date :** 19/01/2024  
**Contact Number :** [Number]  
**Designation :** SUPERINTENDENT  
**Hospital Seal :** DISTRICT HOSPITAL, WADAKKANCHERY, KUMARANELLUR PO THRISUR-680 590

**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT  
MAINTENANCE PROGRAM (BEMP)**

Recommendations for Beyond Economic Repair (BER)

PROFORMA

Sl No	Particulars	Remarks
1	Name of District	THRISSUR
2	Name of Hospital	DH WADAKKANCHERY
3	Name of Equipment with Make, Model and Serial Number	Dialyser reprocessor unit Cura DCS Plus 065/02/18
4	Equipment ID/ Barcode	0812018
5	Date of purchase/ Year of manufacture/Installation Date	10/7/18
6	Warranty details	No Warranty
7	*AMC/ CAMC Period agreed at the time of purchase	No AMC/CMC
8	Date of breakdown/ Date of registration of complaint	17/1/23
9	Action taken	Checked the equipment and found Main board and solenoid valve set defective
10	Present status of the equipment (Fully damaged / partially damaged)	Fully Damaged
11	Details of service required (Recommendations for repair)	Not recommending
12	Cost of spares (specify parts and cost)	NA
13	Asset Value	2,62,500/-


14	Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	NA
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX.	Cyrix Service report attached
16	Reasons for recommending the equipment as BER	Checked equipment and found that Main board Solenoid valve set. Defective.. Equipment installed on 10/07/2018 and aged up to 6 years. Enquired for spare and OEM not found. As per tender clause 5.3.14.2 we are highlighting this case as a exceptional case and recommend equipment for condemnation.
17	Signature of CYRIX Authority	Ramdas K 

\*Not mandatory #Based on the period of life and value as per the BER guidelines


\* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:

Care service not available from company

  
Signature of JC BM (NHM)

Date 31.01.2024

  
Signature of Superintendent/Medical Officer (IC)



SUPERINTENDENT  
DISTRICT HOSPITAL, WADAKKANCHERY  
KUMARANELLUR PO THRISSUR-680 530

10.7.18 Invoice No IN/03/0045/18-19 dt 28.6.18 ref  
 IN/03 Curves Healthcare Pvt. Ltd No. DP 531  
 3rd Main Road Sitalco, Indanur, TUM.

	Qty	Amount
1 DCS plus Dialyzer Repra Cusong Systems.	No	262500.00
1 GST output		31500.00
(Prices Temp labels Ninely four thousand only)		294000.00
		<u>294000</u>

294000

5/2/21 led to Dialyzer Unit

... items correctly  
 No. ....  
 Register Vol. ....  
 Date: .....

cl. to Sub Street  
 Post: 25  
 Diabas

*[Signature]*  
 21/01/2024



SUPERINTENDENT  
 DISTRICT HOSPITAL, SITALCO, TUMKUR  
 KARNATAKA

02/12/17

**CYRIX** MEMORIAL EQUIPMENT MAINTENANCE  
PROGRAM GROUP

TOLL FREE NO: 1800 - 425-7900

BAR CODE - (8004-999818228)  
0812018



**CURA**

WWW.CURA.IN  
Growth. Innovation. Search.

TOLL  
FREE **1800-419-0247**

*Handwritten:* 10/1/2014

*Handwritten:* SU...  
1570...  
1021E-220





MODEL : DCS-PLUS  
Serial no : 065/02/18  
Input : 220± 10V~, 50Hz  
Rated load : 2.4A



Tuscano Equipments Private Limited  
New no: 121, Old no.61,  
First Floor, "Jammi Buildings"  
Royapettah High Road, Mylapore  
Chennai - 600 004, India.

CE



Wellkang Ltd( [www.CE-marking.eu](http://www.CE-marking.eu) )  
29 Harley st, London W1G 9QR, UK  
Ph: +44(20)32876300,  
Fax: +44(20)76811874  
Email: [AuthRep@CE-marking.eu](mailto:AuthRep@CE-marking.eu)



*Handwritten signature and date 3/10/1*

SUPERINTENDENT  
DISTRICT HOSPITAL WADAKKANCHERY  
KUMBAKONAM DISTRICT TAMIL NADU 600 531

CURA



DCS plus

Control panel with a numeric keypad (0-9) and function buttons: START (red arrow), STOP (red arrow), and a red emergency stop button.

**WARNING!**

1. Please read the manual before using the machine. Do not use the machine if you are not familiar with the operation. Do not touch the machine when it is running.

2. Do not touch the machine when it is running. Do not touch the machine when it is running. Do not touch the machine when it is running.

3. Do not touch the machine when it is running. Do not touch the machine when it is running. Do not touch the machine when it is running.

4. Do not touch the machine when it is running. Do not touch the machine when it is running. Do not touch the machine when it is running.



Handwritten signature and date: 3/10/11/024

Official stamp: SUPERINTENDENT  
HOSPITAL