


458

**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAM (BEMP)**

Recommendations for Beyond Economic Repair (BER)

PROFORMA

Sl. No	Particulars	Details
1	Name of District	PALAKKAD
2	Name of Hospital	THIRU ALATHUR
3	Name of Equipment with Make, Model and Serial Number	DIGITAL BP APPARATUS MAKE: LIFECARE MODEL: NIL S.NO: NIL
4	Equipment ID & Barcode	75810, 0920915
5	Date of purchase / Year of manufacture / Installation Date	26-11-2018
6	Warranty details (Yes/No)	NO WARRANTY
7	*AMC/ CAMC Period agreed at the time of purchase	NO AMC / NO CAMC
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	20/10/2023, TOLL FREE
9	Action taken	Checked and found that Display main board, NBP cuff and 4 numbers of batteries faulty
10	Present status of the equipment (Fully damaged / partially damaged)	Fully damaged
11	Recommendations for repair (required service details)	Not recommending
12	Cost of spares (specify parts and cost)	NP


13	Asset Value	2200
14	# Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	NP
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	CYRIX SERVICED REPORT
16	Reasons for recommending the equipment as BER	Checked and found that LED display, main board, BP bulb, BP cuff, BP bladder, BP valve and three numbers of batteries faulty. Machine installed on 26/11/2018 and outlive five plus years. Spare quote not submitted. Spare are not available in market and since recommending for condemnation.
17	Name & Signature of CYRIX Authority	 DIVINIK

*Not mandatory #Based on the period of life and value as per the BER guidelines


* Attach Photograph


Remarks and Recommendations of Junior Consultant (Biomedical) NHM:

Main board, LED display, BP bulb & cuff damaged. Hence recommended for BER.


 Consultant (Biomedical Engg)
 Palakkad
 Signature of J.C. BM (NHM)

Date




SUPERINTENDENT
T.H.O. HOSPITAL
ALATHUR
 Superintendent / Medical Officer (i/c)



**BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAMME
UNDER
NATIONAL HEALTH MISSION**



SERVICE PROVIDER
Tender No. WO-37/2021-2022/698

CYRIX
HEALTHCARE PVT LTD

No. : 164118

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

Service Report 30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

Health Facility : THOH ALATHUR	Call Registration Date : 20/10/23
Address : ALATHUR	Caller ID : 75810
PALAKKAD	Date of Visit : 22/10/23
Ph : 9446139336	Asset No. : 0920915
	EQPT Name : DIGITAL BP APPARATUS
	Manufacture : LIFECARE Model : NIL
	S. No. : NIL Dept. : LABOUR ROOM

Service Classification : Breakdown Call PMS Calibration Cust.Training

Problem Identified : Display defective

Action Taken : Checked the machine. Found that Display, Main board, Air pump, NIBP cuff and 4 numbers of batteries defective. Need to replace these spare for further checking working condition of equipment

Completed Date : 31/10/23 Time : 1:30 PM Spare Required

Spare Replaced Requested


Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrix Engineer	Date	Start Time	End Time
DIVIN'S	31/10/23	12:30 PM	1:30 PM

Customer Remark Completed Pending

Service Engineer Name : DIVIN'S
Signature :
Date : 31/10/23
Contact Number : 7593847155

Customer Name : B.I.V.N.R
Signature :
Date : 31/10/23
Contact Number : 9446139336
Designation : SNO
Hospital Seal :



Superintendent
T.H.Q. HOSPITAL
ALATHUR

B. P. APPARATUS - DIGITAL

		Rd	Issued	Balance	Main stock Balance	Sifa
8/3/21	Rd From the store (80093940	2	✓	2 ✓		Store
9/3/21	Issued Barcode 80093939		2	NIL	2	
17/4/21	B.F From Dage No: 132 Old Back B1	3	3	NIL	5	Store
	(Rd 23/1/2018 - 2 Pins Polite Desiyathu Dem Rs 2000 per 1 Nos, 26/11/2018 - 1					
27/1/23	Issued from T. BOLATION WARD C 1187.	1			6	
9/11/23	Rd from store	1		1	7	BM
10/11/23	Issued for use		1	NIL	7	BM



Jayanthi K Jayath
 Sr SNO BJNR



Jayanthi
 For SNO Biji-N-N