


**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT  
MAINTENANCE PROGRAM (BEMP)**

**Recommendations for Beyond Economic Repair (BER)**

**PROFORMA**

Sl. No	Particulars	Details
1	Name of District	ERNAKULAM
2	Name of Hospital	THQH FORTKOCHI
3	Name of Equipment with Make, Model and Serial Number	SUCTION APPARATUS SURGIVAC C121019
4	Equipment ID & Barcode	63951, 0420734
5	Date of purchase / Year of manufacture / Installation Date	23/4/2005
6	Warranty details (Yes/No)	NO
7	*AMC/ CAMC Period agreed at the time of purchase	9/9
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	9/8/2008
9	Action taken	Checked MOTOR WORTHY CONDITION
10	Present status of the equipment (Fully damaged / partially damaged)	Fully damaged
11	Recommendations for repair (required service details)	Found Suction Appatus MOTOR WAS FULLY RUSTED AND IT CANNOT BEEN REPAIRED
12	Cost of spares (specify parts and cost)	nil


13	Asset Value	8726
14	# Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	Nil
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	CYRIX Service Report
16	Reasons for recommending the equipment as BER	Motor was fully rusted and sector appus was also fully damaged, it cannot be serviced
17	Name & Signature of CYRIX Authority	hiron joshi 


\*Not mandatory #Based on the period of life and value as per the BER guidelines

\* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:


Motor rusted. Completely damaged.  
Submitted BER after 30 days.

  
 CONSULTANT  
 Bio-Medical Engineer  
 National Health Mission (NHM)  
 Arogyakeralam, Mysuru  
 2/12/23  
 Signature of JC BM (NHM)

09  
 Nani P. Jose 

Date

Seal 

  
 Signature of  
 Superintendent / Medical Officer (i/c)

**BIO MEDICAL MAINTENANCE PROGRAMME UNDER NATIONAL HEALTH MISSION**



SERVICE PROVIDER  
Tender No. WO-37/2021-2022/698

**CYRIX**  
HEALTHCARE PVT LTD

No. : 153412

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

**Service Report**

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala  
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

Health Facility : THQH Address : FORTICODU ERNALANDU Ph : 9946755839	Call Registration Date : 19/8/2023 Caller ID : 63951 Date of Visit : 19/8/2023 Asset No. : 02200934 EQPT Name : Suction apphs Manufacture SUPREMAC Model : 719 S. No. Nil Dept. OT
---	--

Service Classification : Breakdown Call  PMS  Calibration  Cust.Training

Problem Identified : Fault Suction apphs not working properly

Action Taken : Checked and found Suction apphs was fully damaged, motor was rusty, it cannot be repaired  
Reorder FOARBE12

Completed  Date : 19/8/23 Time : 10:00am Spare Required

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrix Engineer	Date	Start Time	End Time
h.yeorge	19/8/23	10:am	11:00am

Customer Remark  Completed  Pending

Service Engineer Name : h.yeorge  
 Signature : *[Signature]*  
 Date : 19/8/23  
 Contact Number : 9847299500

Customer Name : Mini P. Jose  
 Signature : *[Signature]*  
 Date : 19/8/23  
 Contact Number : 9946755839  
 Designation : SMO  
 Hospital Seal :

CYRIX BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAM (BEMP)

TOLL FREE NO: 1800 • 425-7669

BAR CODE (800) 590615225  
0720734




14 14 JCNK P 54 10 12 110706

# RATUS TYPE 1

- FIRST IS THE COLLECTION JAR AND SECOND IS FOR OVERFLOW. FLUID SHOULD NOT ENTER THE PUMP/FILTER
- USE THREE-PIN PLUG. ALL ELECTRICAL FITTINGS MUST BE WELL EARTHED
- USE STABILIZER IN CASE OF VOLTAGE FLUCTUATION
- INSTRUCTION MANUAL TO BE FOLLOWED TO MAINTAIN WARRANTY
- SUPPLY AC ONLY 50 HZ
- VOLTAGE 230V • CURRENT 2.6 AMP

C-121-019



BE CAREFUL DRAINAGE OF BODY CAVITIES  
THIS CABINET FOR SAFETY REASONS.



STOCK BOOK OF

Name of Article. SECTION APPROPRIUS

DATE	No. & DATE OF VOUCHER OR INVOICE	BROUGHT FORWARD	RECEIPT	CONSUMED	BALANCE AFTER EACH TRANSACTION	INITIALS OF RECEIVING	REMARK
	5/04/2007	From old Register Vol II A, P.No: 192	1		1		
	18/10/2003	Ret from store (Cao: 123-003)	1		2		
	28/04/2005	Ret from store (Truly model 5.5 Top) (MS P.No: 114)	1		3		
	1/07/2005	Transf to Pre-OP. ward		1	2		
	Old SB 6 pr.wd)				3		
	16/03/2022	To RO from CH ward (Stock page 35) portable	1				