


**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT  
MAINTENANCE PROGRAM (BEMP)**

**Recommendations for Beyond Economic Repair (BER)**

**PROFORMA**

<b>Sl. No</b>	<b>Particulars</b>	<b>Details</b>
1	Name of District	KANNUR
2	Name of Hospital	PHC EZHOME
3	Name of Equipment with Make, Model and Serial Number	SEMI AUTOBIOCHEMISTRY ANLYZER Make: ROBONIK Model: PRIETEST TOUCH SN: AT3700614RBK
4	Equipment ID/ Barcode	1342587 #92848
5	Date of purchase/ Year of manufacture/Installation Date	16-09-2014
6	Warranty details (Yes/No)	No warranty
7	*AMC/ CAMC Period agreed at the time of purchase	No AMC/CAMC
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	Toll free - 26-01-2024
9	Action taken	Checked the machine and found that Low gain showing for all filters found problem with mainboard and filter wheel assembly , enquired OEM for quotation
10	Present status of the equipment (Fully damaged / partially damaged)	Fully damaged
11	Recommendations for repair (required service details)	No recommendations
12	Cost of spares (specify parts and cost)	N/A




13	Asset Value	64760/-
14	# Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	N/A
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	Cyrix Service report & OEM End of life letter attached
16	Reasons for recommending the equipment as BER	Checked the machine and found that Low gain showing for all filters found problem with mainboard and filter wheel assembly, enquired OEM for quotation. Equipment installed on 16-09-2014 aged up to 8 years and 8 months, oem informed that this machine is obsolete model, end of life letter issued. As per tender clause 5.3.14.2 we recommending the equipment for condemnation.
17	Name & Signature of CYRIX Authority	SARANG K M 

\*Not mandatory      #Based on the period of life and value as per the BER guidelines

\* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:


checked & identified that machine model obsolete model. so recommended for BER as per BER protocol.


  
15/2/2024  
Signature of JC BM (NHM)

CONSULTANT BIOMEDICAL  
NATIONAL HEALTH MISSION  
KANNUR - 670002

Namitha.V.K  
~~Namitha~~  
Lab technician grade II  
Block Phc ezhome

Date: 27/02/2024



  
KS  
Dr. Khadees Beebi  
27/2/24.

Signature of  
Superintendent/Medical Officer (i/c)

ബി.എ.ടി. മൂന്നാം ഫിസിക  
എസ്.എം. കെ.എസ്. 670 334  
ഫോൺ: 0497 2815163





# BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME UNDER NATIONAL HEALTH MISSION



**SERVICE PROVIDER**  
Tender No. WO-37/2021-2022/698

## CYRIX

HEALTHCARE PVT LTD

No. : 174311

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

**Service Report** 30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala  
Ph : 98472 99500 Website : www.cyrrix.com | E-mail : bemp.kl@cyrix.in

<p>Health Facility ..... <u>PHC Ezhome</u></p> <p>Address ..... <u>Kannur</u></p> <p>Ph : ..... <u>8891509784</u></p>	<p>Call Registration Date : <u>26-01-2024</u></p> <p>Caller ID : <u>92848</u></p> <p>Date of Visit : <u>27-01-2024</u></p> <p>Asset No. : <u>1342587</u></p> <p>EQPT Name : <u>Semi Auto Biochemistry Analyzer</u></p> <p>Manufacture <u>Robonik</u> Model : <u>Prietest Touch</u></p> <p>S. No. <u>AT370064RFX</u> Dept. <u>Laboratory</u></p>
---	---

**Service Classification :** Breakdown Call  PMS  Calibration  Cust.Training

**Problem Identified :** Found low gain showing for all filters.

**Action Taken :** Checked the machine and found that low gain showing for all filters found problem with mainboard and filter wheel assembly, enquired OEM for quotation.

Completed  Date : 27-01-2024 Time : 11:00AM Spare Required

Spare Replaced  Requested

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrix Engineer	Date	Start Time	End Time
<u>Nitheesh.P</u>	<u>27-01-2024</u>	<u>10:00AM</u>	<u>11:00AM</u>

**Customer Remark** Completed  Pending

*Dr. What's app's Rehi 27/2/24*

<p><b>Service Engineer Name :</b> <u>Nitheesh.P</u></p> <p><b>Signature :</b> </p> <p><b>Date :</b> <u>27-01-2024</u></p> <p><b>Contact Number :</b> <u>8301871591</u></p>	<p><b>Customer Name :</b> <u>Alamitha</u></p> <p><b>Signature :</b> </p> <p><b>Date :</b> <u>27/01/24</u></p> <p><b>Contact Number :</b> <u>9562728983</u></p> <p><b>Designation :</b> <u>Lab technician grade II</u></p> <p><b>Hospital Seal :</b></p>
--	---




ഹെൽത്ത് സർവീസ് ഓഫീസർ  
 ഹെൽത്ത് സർവീസ് ഓഫീസർ  
 ഹെൽത്ത് സർവീസ് ഓഫീസർ  
 ഹെൽത്ത് സർവീസ് ഓഫീസർ  
 ഹെൽത്ത് സർവീസ് ഓഫീസർ



**Product BIOCHEMISTRY ANALYSER** **IVD**

Model **prietest TOUCH** Voltage **115V-230V/ 60Hz-50Hz**

Sr. No. **AT3700614RBK** Power **75VA** **CE**



*K...*  
*22/2/24*

കുമാരൻ  
 സഹായക  
 സേവനം  
 ഫോൺ: 0471 2815163

Date: 30<sup>th</sup> January 2023

Dear Valued Customer,



Thank you for choosing **Robonik (India) Pvt Ltd**. We appreciate the opportunity to support our business. This letter is to inform you regarding the EOL (End of Life) of the instrument Semi Automated Biochemistry Analyser (Model: Prietest Touch) having **Serial No. AT3700614RBK**.

As required, we stated herewith that the said instrument Semi Automated Biochemistry Analyser (Model: Prietest Touch) having **Serial No. AT3700614RBK** as **End of Life (EOL) Notified**. (Duration from installation is 10 Years).

Thank you in advance for your understanding and your cooperation. If you have any questions, please contact our Customer Care Department at Toll Free No. 18005727977.

Yours faithfully,

For **ROBONIK INDIA PVT LTD**



**AUTHORIZED SIGNATORY**



# Stock Register of

Name of Drug: **EQUIPMENT SEMI AUTO ANALYSER**

Date	Particulars of Receipt or Issue	Receipt	Issue	Balance
1	2	3	4	5
21.05.15	INNO: 1/015-16 CHC Perigone MODEL Pyetest To web	1 NO		1 NO
16.06.15	Issued to lab		1 NO	NIC



# Medicines

212

Unit	Rate		Cost	Batch No.	Name of Manufacturer	Date of Expiry	Initial
	Rs.	Ps.					
6	7		8	9	10	11	12
1			64760		Robonic		<del>Handwritten initials</del>

