


**REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT  
MAINTENANCE PROGRAM (BEMP)**

**Recommendations for Beyond Economic Repair (BER)**

**PROFORMA**

<b>Sl. No</b>	<b>Particulars</b>	<b>Details</b>
1	Name of District	KANNUR
2	Name of Hospital	FHC ALAKODE(T)
3	Name of Equipment with Make, Model and Serial Number	WEIGHING MACHINE Make :HEALCURE Model:N/A SN:N/A
4	Equipment ID/ Barcode	1340813 #104358
5	Date of purchase/ Year of manufacture/Installation Date	27-05-2018
6	Warranty details (Yes/No)	No warranty
7	*AMC/ CAMC Period agreed at the time of purchase	No AMC/CAMC
8	Date of breakdown (Date of registration of complaint through email/ Toll free)	Toll free- 02-04-2024
9	Action taken	Checked the machine and found that Machine is not switching on, Found main board, Load cells damaged. Need to replace the these spares for further checking .
10	Present status of the equipment (Fully damaged / partially damaged)	Fully damaged
11	Recommendations for repair (required service details)	No recommendations
12	Cost of spares (specify parts and cost)	NA

13	Asset Value	850 /-
14	* Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	NA
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not)	Service report attached
16	Reasons for recommending the equipment as BER	Checked the machine and found that machine is not switching ON. Found that main board and load cells damaged. Equipment installed on 27.05.2018 aged up to 5 years and 11 months. Quotation not submitted since Spare parts are not available on market. So we recommend the equipment for condemnation.
17	Name & Signature of CYRIX Authority	AVINASH T 


\*Not mandatory

#Based on the period of life and value as per the BER guidelines


\* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:

checked & identified that machine not switching on. its main board & load cells are damaged. aged up to 5 years & above. so recommended for BER as per BER protocol.

  
17/04/24.

Signature of JC BM (NHM)

Srjg. C. U   
Nsg officer

Date: 18/4/2024



  
Signature of  
Superintendent Medical Officer (S/O)

(Kannada text)

(Kannada text)





**BIOMEDICAL EQUIPMENT  
MAINTENANCE PROGRAMME  
UNDER  
NATIONAL HEALTH MISSION**



**SERVICE PROVIDER**  
Tender No. WO-37/2021-2022/698

**CYRIX**  
HEALTHCARE PVT LTD

No. : 185636

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

**Service Report**

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala  
Ph : 98472 99500 Website : www.cyrrix.com | E-mail : bomp.kl@cyrrix.in

Health Facility <u>FHC ALAKODE (T)</u>	Call Registration Date : <u>2/4/24</u>
Address <u>THERTHALLY</u>	Caller ID : <u>104398</u>
<u>KANNUR</u>	Date of Visit : <u>3/4/24</u>
Ph : <u>9745939544</u>	Asset No. : <u>1340813</u>
	EQPT Name : <u>weighing machine</u>
	Manufacture <u>healux</u> Model : <u>na</u>
	S. No. <u>na</u> Dept. <u>OP</u>

**Service Classification :** Breakdown Call  PMS  Calibration  Cust.Training

**Problem Identified :** Machine is not switching ON

**Action Taken :** checked the machine and found that machine is not switching ON. found main board lead cells damaged. Need to replace the these spares for further checking

Completed  Date : 3/4/24 Time : 11:10am Spare Required

Spare Replaced  Requested

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrrix Engineer	Date	Start Time	End Time
<u>Aswin</u>	<u>3/4/24</u>	<u>10:30am</u>	<u>11:10 am</u>

**Customer Remark** Completed  Pending

<p><b>Service Engineer Name :</b> <u>Aswin</u></p> <p><b>Signature :</b> <u>[Signature]</u></p> <p><b>Date :</b> <u>3/4/24</u></p> <p><b>Contact Number :</b> <u>8301012759</u></p>		<p><b>Customer Name :</b> <u>Suje - C.U</u></p> <p><b>Signature :</b> <u>[Signature]</u></p> <p><b>Date :</b> <u>3/4/24</u></p> <p><b>Contact Number :</b> <u>9745939544</u></p> <p><b>Designation :</b> <u>NSG officer</u></p> <p><b>Hospital Seal :</b> <u>[Signature]</u></p>
---	---	--

WEIGHING MACHINE

	R	J	39	890	890
25/1/22 in (new account book volume I) Page No. 154 B/F → (Digit 2 + ordinary records)			3	2	2
4-1-2019 Condensed (No ordinary) as per order No 365/18 dated 4-1-22 of M.O of FHC thettdly.		1	2	2	(Ord type)
21/12/22 Reviewed from the store (ordinary type) RS -	1		3	2	



Kumar EV

സെക്ഷൻ ഓഫീസർ  
കെ.പി.എസ്. ഓഫീസ് കെ.എസ്.എസ്.



CYRIX  
TOLL FREE NO: 1800 - 423-7889  
MAX CODE: (800-800-1522)  
1340813



Max. 150kg d=100g





ИЗДЕЛІЕ СЕРИЯЛІЗІРОВАНО  
 CE

MESL KTRPL

00441 800 7479 - 8710016



**MES**

Approved	01.07.2016
By	...
For	...
At	...
Location	...

Manufacturing & Service  
 MESL KTRPL

**CYRIX**

СЕРІЯЛІЗІРОВАНО  
 ЦИФРОВИЙ ІДЕНТИФІКАЦІЙНИЙ  
 СИСТЕМА

Barcode



Копія документа, наданого на вимогу  
 Інформаційно-технологічного департаменту  
 - 67 N 5711-16-11