


Ob
8/10/2

REPAIR OF BIOMEDICAL EQUIPMENTS UNDER BIOMEDICAL EQUIPMENT
MAINTENANCE PROGRAM (BEMP)

Recommendations for Beyond Economic Repair (BER)

PROFORMA

Sl. No	Particulars	Details
1	Name of District	THRISSUR
2	Name of Hospital	DH WADAKKANCHERRY
3	Name of Equipment with Make, Model and Serial Number	Dialyzer Reprocessing Unit MAKE: Cura MODEL: DCS plus SN: 065/02/18
4	Equipment ID/ Barcode	0812018 # 91214
5	Date of purchase/ Year of manufacture/Installation Date	10-07-2018 (stock book attached for reference)
6	Warranty details	No Warranty
7	*AMC/ CAMC Period agreed at the time of purchase	No AMC/CAMC
8	Date of breakdown/ Date of registration of complaint	17-01-2024
9	Action taken	Checked the machine and found that Main board, solenoid and display valve set are defective. need to replace the solenoid valve set, display and main board for further checking.
10	Present status of the equipment (Fully damaged / partially damaged)	Fully Damaged
11	Details of service required (Recommendations for repair)	Not recommending
12	Cost of spares (specify parts and cost)	NA


13	Asset Value	2,62,500/-
14	# Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value	NA
15	Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX.	Cyrix Service report attached
16	Reasons for recommending the equipment as BER	regarding the inspection of the equipment, dialyzer reprocessing unit. Upon examination, it has been discovered that critical spare parts, including the Main board, Solenoid valve and display are severely damaged. Equipment installed on 10/07/2018 and aged up to 5+year In our effort to source replacement parts, we inquired with the original manufacturer, Cura Healthcare. However, it has come to our attention that Cura Healthcare has ceased its operations and has permanently closed down. Consequently, there are no available spares from the manufacturer. Considering the circumstances and in adherence to the tender clause 5.3.14.2, which stipulates that manufacturer not existing, the equipment should be evaluated for condemnation, we recommend that these equipment is marked for condemnation.
17	Signature of CYRIX Authority	BIJO T JOY 

*Not mandatory #Based on the period of life and value as per the BER guidelines

* Attach Photograph

Remarks and Recommendations of Junior Consultant (Biomedical) NHM:
Company closed and no spare support

[Signature]
 Signature of JC BM (NHM)
 Jr. Consultant (Biomedical)
 NHM, Thrissur

Superintendent/Medical Officer (i/c)  Signature of seal Date *[Signature]*

SUPERINTENDENT
 DISTRICT HOSPITAL, WADAKKANCHERY
 KUMARANELLUR PO. THRISSUR - 680 590

10.7.18 Invoice No 1N/03/0045/18-19 dt 28.6.18 ref
 1N/03 Curva Healthcare Pvt. Ltd No. DP 561
 3rd Main Road Sidco, Indaravadi, T.M.

Qty	Amount
1 DCS plus Dialyzer Repra cussing system	262500.00
1 GST output	31500.00
(Proposed Temp labels Novelty four thousand only)	294000.00
	<u>294000</u>

294000

5/2/21 led to ... items correctly
 Unit

Signature: _____
 Name: _____
 Designation: _____
 Date: _____

cl to substance
 pmt: 25
 Dinkus

[Signature]
 2/10/2024



SUPERINTENDENT
 DISTRICT HOSPITAL WADAKKANCHEERY
 KULMRANELUR PO THIRISSUR. 688 530



BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME UNDER NATIONAL HEALTH MISSION



SERVICE PROVIDER
Tender No. WO-37/2021-2022/698

CYRIX

HEALTHCARE PVT LTD

No. : 178906

ISO 13485 : 2012 & ISO 9001-2008 CERTIFIED COMPANY | AERB Approved Service Agency

Service Report

30/64 1 B, Petta Junction, Poonithura, Kochi - 682 038, Kerala
Ph : 98472 99500 Website : www.cyrix.com | E-mail : bemp.kl@cyrix.in

<p>Health Facility <u>DH Wadakkanchery</u></p> <p>Address <u>Thrissur</u></p> <p>Ph : <u>7560935204</u></p>	<p>Call Registration Date : <u>17/01/2024</u></p> <p>Caller ID : <u>91214</u></p> <p>Date of Visit : <u>19/01/2024</u></p> <p>Asset No. : <u>0812018</u></p> <p>EQPT Name : <u>DIALYSER REPROCESSOR UNIT</u></p> <p>Manufacture <u>CURA</u> Model : <u>DOS PLUS</u></p> <p>S. No. <u>06510218</u> Dept. <u>DIALYSIS</u></p>
---	---

Service Classification : Breakdown Call PMS Calibration Cust.Training

Problem Identified : Not working

Action Taken : checked the Equipment found that Main board and Solenoid valve set defective.

Completed Date : Time : Spare Required

Spare Replaced Requested

Description	Qty.	Part Number	PR Number
1.			
2.			
3.			

Cyrix Engineer	Date	Start Time	End Time

Customer Remark Completed Pending

Service Engineer Name : Tejas

Signature : [Signature]

Date : 19/01/2024

Contact Number : 9848098422

Customer Name : [Signature]

Signature : [Signature]

Date : [Date]

Contact Number : [Number]

Designation : SUPERINTENDENT

Hospital Seal :

DISTRICT HOSPITAL WADAKKANCHERY
KUMARANELLUR PO THRISSUR-680 590



0815512

CYRIX
CARE PVT LTD
S.N:065
Reprocessing Unit
2023
2024
Certified Co. Ph: 9744555073

CYRIX BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME (BEMP)
TOLL FREE NO: 1800 - 425-7669
BAR CODE - (8004-890615225) 0812018




CURA
www.cura.in
Growth. Innovation. Reach

ADK
31/01/2024

TOLL FREE **19-0247**

SUPPORT CENTER
DISTRICT HOSPITAL WADAKKANCHERY
KUMARANELLUR PO THIRUVANANTHAPURAM





MODEL : DCS-PLUS
 Serial no : 065/02/18
 Input : 220±10V~, 50Hz
 Rated load : 2.4A



Tusciano Equipments Private Limited
 New no: 121, Old no.61,
 First Floor, "Jammi Buildings"
 Royapettah High Road, Mylapore
 Chennai - 600 004, India.



Wellkang Ltd(www. CE-marking.eu)
 28 Harley st, London W1G 9QR, UK
 Ph: +44(20)32876300,
 Fax: +44(20)76811874
 Email: AuthRep@CE-marking.eu




Handwritten signature
 3/10/1

SUPERINTENDENT
 DISTRICT HOSPITAL, WADAKKANCHERY
 KUMARANELLUR PO THRISSUR, 680 500

CURA

DCSplus



START ↑

STOP ↓

ENTER

1 2 3

4 5 6

7 8 9

0

BLOOD PORT IN

WARNINGS!

Read INSTRUCTION MANUAL prior to using this machine (Dialyser Cleaning System). Failure to follow the instructions for use may result in injury.

If a Dialyser is processed on DCSplus with a test dialyser, care must be taken to use CHEMICAL DIALYSATE. For safety, a special test dialyser must be used.

DCSplus, Dialyser Cleaning System is not to be used for the purpose of dialysis. Always marking that a dialyser is used for dialysis. Do not use DCSplus for dialysis. Always marking that a dialyser is used for dialysis. Do not use DCSplus for dialysis.

1) Do not use DCSplus for dialysis.

2) Do not use DCSplus for dialysis.

3) Do not use DCSplus for dialysis.

4) Do not use DCSplus for dialysis.

5) Do not use DCSplus for dialysis.

6) Do not use DCSplus for dialysis.

7) Do not use DCSplus for dialysis.

8) Do not use DCSplus for dialysis.

9) Do not use DCSplus for dialysis.

10) Do not use DCSplus for dialysis.

11) Do not use DCSplus for dialysis.

12) Do not use DCSplus for dialysis.

13) Do not use DCSplus for dialysis.

14) Do not use DCSplus for dialysis.

15) Do not use DCSplus for dialysis.

16) Do not use DCSplus for dialysis.

17) Do not use DCSplus for dialysis.

18) Do not use DCSplus for dialysis.

19) Do not use DCSplus for dialysis.

20) Do not use DCSplus for dialysis.

BLOOD PORT OUT



Handwritten signature
31/01/2024

SUPERINTENDENT
DISTRICT HOSPITAL WADAKKANCHERY