

PROFORMA
RECOMMENDATION FOR BEYOND ECONOMIC REPAIR (RBER)
BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAM (BEMP)

Name of Hospital: PHC KOLACHERY	Name of District: KANNUR
Name of Equipment: BIOCHEMISTRY ANALYSER	Equipment ID & Barcode: 136301-1342439
Make: ROBONIK	Model: PRIETEST TOUCH
Serial Number: AT0820815RBK	Warranty details: No Warranty
Date of purchase/ Year of manufacture/ Installation Date: 01-03-2016	Present status of the equipment: Fully damaged
Date of breakdown: 23.09.2024 Date of registration of complaint through email/Toll free): Toll free	*AMC/ CAMC Period agreed at the time of purchase: No AMC/CAMC

Action taken: Checked the Equipment Found Equipment Is not Switching ON. Found That primary Board and Display Defective. Enquired For quotation.

Recommendations for repair (required service details): No Recommendations

Cost of spares (specify parts and cost): NA

#Percentage value of the cost of spares with respect to Cost of Purchase/ Asset Value: NA	Asset Value: 64,760
---	---------------------

Abstract of Service Report provided by the OEM/ Authorized Service Provider/ CYRIX (Attached or Not): Cyrix Service Report & OEM end of life Letter attached.

Reasons for recommending the equipment as BER Checked the Equipment Found Equipment Is not Switching ON. Found That primary Board and Display Defective. Enquired For quotation. Equipment Installed On 01.03.2016, aged up to 8 years 6 months. OEM Informed that this machine is obsolete model, end of life letter issued. As per Tender Clause 5.3.14.2 we recommend the Machine for condemnation.

Avinash T Amal
 Name & Signature of CYRIX Authority with date 30/09/24

Remarks of Junior Consultant (Biomedical) NHM: checked & identified that the machine not switching ON due to Primary board & Display board are defective. Aged up to 8 years & above. The OEM informed that machine is obsolete model. As per BER Protocol recommended for BER.

Recommended for BER (Yes/ No):	<i>Rathin</i> 30/09/24
Date: 30/09/24.	Signature of JC BM. (NHM)

Date: 30/9/24	Seal	Signature of Superintendent / Medical Officer (i/c)
---------------	------	---

*Not mandatory

#Based on the period of life cycle as per the BER guidelines

Attach Photograph



CONSULTANT BIOMEDICAL
 NATIONAL HEALTH MISSION
 കണ്ണൂർ ജില്ലാ ആരോഗ്യ മന്ത്രിയുടെ കാര്യാലയം
 ഫോൺ - 670604

Date: 27.09.2024

Dear Valued Customer,

Thank you for choosing **Robonik (India) Pvt. Ltd.** We appreciate the opportunity to support our business. This letter is to inform you regarding the EOL (End of Life) of the instrument Semi Automated Biochemistry Analyser (Model: Prietest Touch) having **Serial No. AT0820815RBK.**

As required, we stated herewith that the said instrument Semi Automated Biochemistry Analyser (Model : Prietest Touch) having **Serial No. AT0820815RBK** as **End of Life (EOL) Notified.** (During from installation is 09 years).

Thank you in advance for your understanding and your cooperation. If you have any questions, please contact our Customer Care Department at Toll Free No. 18005727977.

Yours faithfully,

For **ROBONIK (INDIA) PVT. LTD.**

A handwritten signature in blue ink is written over a circular blue stamp. The stamp contains the text 'ROBONIK (INDIA) PVT. LTD.' around its perimeter.

Authorised Signatory

+919833696980

EOL Letter = Prietest Touch - AT2311214RBK.pdf

1 page · 89 kB · PDF

11:13 am


Thank uu 11:13 am



11:13 am


Friday

Forwarded



12:58 pm

Forwarded



EOL Letter = Prietest Touch - AT0820815RBK.pdf

1 page · 89 kB · PDF

3:48 pm

Thank you 3:49 pm

You're welcome 3:50 pm

Message deleted for me

UNDO

Message



CYRIX BIOMEDICAL EQUIPMENT MAINTENANCE PROGRAMME (BEMP)

TOLL FREE NO: 1800 - 425-7600


BAR CODE - (8004-800815225)
1342439

CYRIX
HEALTHCARE EQUIPMENT

CN: DMKQ S.N.
PRO: Biochemistry Daily 100
ID: 1000000
CALON: 10-12-20
CALHE: 10-12-20



CYRIX HEALTH CARE PROGRAMME (BEMP) SERVICE PROVIDED THROUGH THE MAINTENANCE PROGRAMME



Product - Biochemistry Analyser IVD
Model: prietest TOUCH Voltage: 115-230VAC, 50-60Hz
SN: AT0820815RBK Power: 100 VA CE

Jinsha T.K

Technician

